







# Quingo™

## Feature Comparison v4

NB. Features/specifications may vary region to region

	PORTABLE		PERFORMANCE		ADVANCED	
	Quingo Air 2	Quingo Flyte & Docking Station	Quingo Classic	Quingo Plus	Quingo Vittess 2	Quingo Toura 2
						
Max Speed Where Permitted	4 mph 6.4 kph	4 mph 6.4 kph	4 mph 6.4 kph	8 mph 12.9 kph	8 mph 12.9 kph	10.5 mph 17 kph
Carrying Capacity	300 lbs 136 kg 21.5 st	350 lbs 159 kg 25 st	350 lbs 159 kg 25 st	350 lbs 159 kg 25 st	400 lbs 182 kg 28.6 st	485 lbs 220 kg 35 st
Maximum Range*	10 miles (opt 18 miles) 16 km (opt 29 km)	23 miles (opt 34 miles) 37 km (opt 54 km)	23 miles 37 km	30 miles 48 km	30 miles 48 km	55 miles 88 km
Turning Radius	104 cm 41 inches	102 cm 40.2 inches	107 cm 42 inches	107 cm 42 inches	113 cm 44.5 inches	132 cm 52 inches
Maximum Kerb Height	5 cm 2"	7.5 cm 3"	7.5 cm 3"	7.5 cm 3"	7.5 cm 3"	10 cm 4"
Maximum Gradient	6 degrees	6 degrees	12 degrees	12 degrees	12 degrees	12 degrees
Battery Type	2 x 15AH (opt 4 x 15AH)	4 x 22AH (opt 6 x 22AH)	2 x 35AH	2 x 50AH	2 x 50AH	2 x 100AH
Motor Size	DC 24 Volt 110 Watt	DC 24 Volt 250 Watt	DC 24 Volt 250 Watt 4 pole	DC 24 Volt 400 Watt 4 pole	DC 24 Volt 400 Watt 4 pole	DC 24 Volt 650 Watt 4 Pole
Off Board Battery Charger	2.0 AH	5.0 AH	5.0 AH	5.0 AH	5.0 AH	8.0 AH (opt 15 AH)
Warranty	Refer to distributor	Refer to distributor	Refer to distributor	Refer to distributor	Refer to distributor	Refer to distributor
Quintell™ 5 Wheel Stability System	✓	✓	✓	✓	✓	✓
Quintell™ Independent Adaptive footplates	✓	✓	✓	✓	✓	✓
Hydraulic Tiller Adjust	✓	-	✓	✓	✓	✓
Quintell™ Kerbmaster Powered Anti-Grounding/ Anti-Tipping	✓	✓	✓	✓	✓	Not Required
Forward Reverse Seat Slider	✓	✓	✓	✓	✓	✓
Quintell™ Smartbumper	-	✓	✓	✓	✓	✓
All Weather Removable Canopy Ready	-	-	✓	✓	✓	✓
Front & Rear Suspension	-	-	✓	✓	✓	✓
Full Front & Rear Lighting Kit	-	-	✓	✓	✓	✓
Easy Grip Manual Disc Brake	-	✓	✓	✓	✓	✓
High Back Seat Support	-	-	-	✓	✓	✓
Self Leveling Tiller Controls	-	-	-	-	✓	✓
EZ View Digital Dash with on screen status/fault indicators	-	✓	-	-	✓	✓
Extended Access Seat & Fold Flat Feature	-	-	-	-	✓	✓
Luxury High Back Seat	-	-	-	-	✓	✓
Ultra Bright LED Low Energy Marker Lamps	-	✓	-	-	✓	✓
12 volt Outlet Socket for Sat Nav / iPod etc	-	-	-	-	✓	✓
Rear View Camera & Heated Handlebar Option available	-	-	-	-	✓	✓
Luxury Extra Wide Comfort Seat with Padded Armrests	-	-	-	-	-	✓
Large Diameter Front Wheel Configuration	-	-	-	-	-	✓
Smooth Drive II Extended Travel Suspension	-	-	-	-	-	✓
Easily portable for transit in Small Cars, Estates & SUV's	✓	✓	-	-	-	-
Remote Control Self-Loading and Docking System	-	✓	-	-	-	-
Range Extender options	✓	✓	-	-	-	-